L'ACCELERATORE DI ELETTRONI MOBILE PER RADIOTERAPIA INTRAOPERATORIA (IORT)  
*THE MOBILE ELECTRON ACCELERATOR FOR INTRA OPERATIVE RADIOTHERAPY (IORT)*

**SORDINA**

tecnologie per la vita / *technology serving life*

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## TECHNICAL IORT PRESENTATION

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*Intraoperative radiotherapy or IORT consists in treating tissue suffering from neoplastic processes, made accessible by surgery, with an electron beam.*

*After an experimental stage, in the sixties last century IORT cleared and set its election criteria and in the mid nineties, thanks to specific equipment enabling validation of new protocols, began to spread outside research environments.*

*IORT has shown an evident capacity for the local control of cancer and can be applied to all solid tumours.*

*The LIAC (Light Intraoperative Accelerator) is a machine created after carefully analysing limits set by what has to be done in a surgical environment; its prerogatives are the excellent possibility of being included in any operating theatre with no need to change existing structures.*

*Its weight, extremely limited (400 kg), enables use on all floor types with no need to reinforce it and thanks to its size (200x80x180 cm) and manoeuvring capacity it can be transported in normal lifts carrying stretchers.*

*The electronic beam energy (10 or 12 MeV) can satisfy almost any kind of therapy of a minimum complexity level to obtain use permits, while the amount of radiation produced means minimum time (40-120 seconds) needed to carry out treatment.*

*The LIAC also has a supply stabilising system providing short and long term stability like conventional accelerators, and its electron beam means execution doses are precise like with completely conventional methods and instruments.*

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*The LIAC was designed to operate in the entire range of intraoperative therapies; it is practical to use in those applications where rapid positioning and considerable radiation volume do not lengthen short operations; however, its energy and agility correspond excellently to the needs of more laborious applications like those on the abdomen and mediastinum.*

*Special care was taken over making therapy easier, where positioning is more delicate; for example the Miles technique.*

*We must point out that the time needed for a radiation therapy is independent of surgery complexity and trained staff do not take longer than 5-7 minutes to supply radiation.*

### SAFETY IN AND OUTSIDE THE OPERATING THEATRE

*A lot of care was taken over machine safety; the position of mobile beam screening is monitored continuously and blocks the machine if it is misaligned. There are far more optical and acoustic signals than the laws in force require.*

*Finally, the LIAC has excellent connectivity as it can be connected to a local network to exchange data; the machine is compatible with all main protocols in use.*

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### EXAMPLE OF HOW A IORT WORKS

**1.** The machine is usually kept in a room next to the operating theatre together with its control cabinet and mobile barriers. About 15 minutes before therapy, the machine is turned on and as soon as the surgeon starts the IORT stage calling the Radiotherapist...

**2.** ...the machine is taken into the theatre; this takes place without turning it off.

**3.** The machine is taken close up to the patient with no need for the latter to be moved; at this stage, the radiotherapist and surgeon decide treatment parameters, fixing the applicator in the surgical cut and documenting the field to be radiated.

**4.** The machine is connected to the applicator.

**5.** While the Health Physicist sets the data for the radiation required.

**6.** The mobile barriers are placed around the table by the radiotherapy technician under the control of the person in charge of radioprotection. The treatment is carried out and lasts between 40 and 120 seconds.

**7.** When treatment is over the machine is taken out of the theatre again.

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<sup>1</sup> On request: 120 mm circular; semicircular and square for composite field; small (2-4-6-8-12-16-20mm) for special applications.

## DIMENSIONI *DIMENSIONS*